



Stable -
Comfortable -
Individual

VACO[®]pedes

VACO[®]pedes - the Safe Solution for your feet

OPED



Stable

- Vacuum cushion and plastic frames effectively stabilize the foot
- Optimal pressure relief in the fore-foot region
- Functional treatment - the talocrural (ankle) joint is freely mobile

Comfortable

- Rocker bottom ensures a good walking gait
- Vacuum cushion provides a comfortable foot bed
- Quick and simple to apply thanks to the Velcro closure system
- Terry cloth covers can be washed and replaced



Individual

- Vacuum cushion adapts to anatomical situation
- Reversible adjustment, anytime, anywhere, without assistance
- Only two sizes; product can be used on either foot: no left, no right version



Indications

- Toe fractures
- Metatarsal fractures
- Hallux valgus/rigidus
- Arthrodeses of the toe joints
- Hammer and claw toes

Available sizes (Note: VACOpedes can be used on either foot)

Standard	Shoe sizes	US men 7-13 / women 9+	EUR 41-46
Small	Shoe sizes	US men 5-7 / women 7-9	EUR 36-40

For more information, please visit our Internet site at www.vacoped.com.

For ordering information, please refer to our price list.



OPED AG • Sumpfstraße 5 • CH-6312 Steinhausen
 Fon +41 (0)41/748 53 35 • Fax +41 (0)41/748 53 40
 E-mail contact@oped.ch • www.oped.ch

OPED Inc. • 47 Loring Drive • Framingham, MA 01702 • USA
 Phone +1/508/879 0762 • Fax +1/508/879 0769 • Toll free +1/866/447 0601
 Email contact@opedinc.com • www.opedinc.com