

More information:  
[www.vacotalus.com](http://www.vacotalus.com)



Keeps you going.

**OPED GmbH**  
Medizinpark 1  
83626 Valley/Oberlaidern  
Germany  
Fon +49(0)8024/60818-210  
Fax +49(0)8024/60818-299  
mail@oped.de  
www.oped.de



WM-INV-T-01-1-EN

**VACO**<sup>®</sup> by OPED  
**talus**

For Stable  
Treatment of the  
Ankle Region



**OPED**  
Keeps you going.

## Advantages

- Optimised combination of effective stabilisation and superb wearing comfort
- Tangible limitation of talar shift through the Talusstrap
- Optimised treatment of swelling since the orthosis offers space to accommodate it
- Excellent grip around the ankle area thanks to VACO12 cushion
- Slim design to allow it to be worn easily inside a shoe

## Indications

- Acute ankle ligament injuries
- After-care for fractures of the ankle joint
- Prevention of ligament lesions on the ankle joint
- Compensation of ankle joint instability

## Product Details

- 1 **Elastic straps**  
for high wearing comfort
- 2 **Limitation of talar shift**  
by the Talusstrap
- 3 **Natural flow of movement of the foot**  
remains possible through plantar/dorsal movement
- 4 **Anatomically pre-moulded outer shell**  
with excellent rigidity for effective stability  
against pronation/supination
- 5 **VACO12 pad**  
for optimised adaptation to individual  
anatomical requirements

